

**Does a QLink work like magnets?**

No. Magnet therapy claims to provide a different benefit and employs different principles.

**How are QLinks different from a magnet therapy product?**

Magnet therapy products claim to relieve and eliminate physical pain on a specific area of the body. QLinks boost your overall mental, emotional and physical energy states at the most fundamental level of life energy.

**What does the term "QLink" mean?**

The term QLink means "quantum-link." Subtle energy science hypothesizes that subtle energies may be a quantum field phenomena, or may influence quantum phenomena.

**What's inside my QLink Pendant?**

The SRT Resonator™ inside your QLink Pendant is made up of three hardware components:

1. The resonating cell. The resonating cell is the size of a watch battery and is made of a unique combination of 99.999% pure crystalline elements.
2. The tuning board. The tuning board is 14 karat gold-plated.
3. The amplifying coil. The amplifying coil is 75 feet of copper wire wound into approximately one inch in diameter.

The main component, the resonating cell, works as a crystalline oscillator to convert and conduct subtle energy into more coherent energy states, ones which are conducive to the healthy energy states of the person. The tuning board reinforces the resonating cell to function in a specific spectrum or energy range. It harmonizes and protects the functioning of the resonating cell.

The amplifying coil shapes and strengthens the size of the energy field that is conducted. It extends the power of the sympathetic resonance into a 3-dimensional field. This entire process strengthens the internal power and resilience of the wearer, allowing him or her to more effectively filter out chaotic energies.

Electrical energy is not required for this SRT Resonator™ since it is powered by the energy of the person wearing it. Because electrical power is not applied, and because of the unique way an SRT Resonator™ works, the cell, tuning board and amplifying coil do not have to be physically wired as they would with a normal Hertzian electronic circuit.

**What's inside my QLink ClearWave or Ally?**

The SRT™ Resonator inside the ClearWave or Ally is a unique combination of components. These are not traditional electrical engineering-type circuit boards, but are designed for their own purpose to convert and conduct subtle energies. These circuit boards use an actively-powered oscillator to extend the influential range of the instrument.

**Where can I get more information about the QLink family of products?**

The QLink product line is manufactured by Clarus Products International (CPI) L.L.C., and can be purchased from CPI or a QLink Independent Distributor by calling 800-246-2765, or by accessing the Clarus web site at [www.clarus.com](http://www.clarus.com).

**How can I learn about QLink business opportunities?**

If you want to learn about QLink business opportunities, call us at 800-246-2765 or send an e-mail to [rw@clarus.com](mailto:rw@clarus.com). A Clarus representative will contact you to discuss your interest and explain this exciting and lucrative opportunity.